## TOP**DENT PERLADIA**

Abrasive rubbers with large filling content of high pure diamonds

Due to the unique bonding technology of rubber and diamonds the PFRI ADIA enables simultaneous removal and polishing of zirconia, pressed ceramics and ceramics in one operation on nano level. Therefore plaque adhesions to ceramic restorations are kept to a minimum.



- Long service life
- Can be used up to the head of the mandrel
- High quality and workmanship
- Efficient removal with simultaneous nano polishing
- · Low heat generation due to high diamond content





Point dark brown Ø 5 mm, L 12 mm REF 20631



dark brown = medium grit

High removal with simultaneous polishing



Point pink Ø 5 mm, L 12 mm REF 20651



Low removal with simultaneous high gloss polishing





